



# **Communication Aids**

### Big Mack – Various colours

COMM 175, COMM 176, COMM 197, COMM 198

Record any single message directly into the Big Mack communicator and press its activation surface for playback up to two minutes in length.



### Little Mack – Various colours

COMM 195, COMM 196

Record any series of messages into the Little Mack. Press its activation surface once to hear the first message. Press it again to hear the next message in the sequence, and so on.



### Go Talk 4, 9 and 20

COMM 080, COMM 081, COMM 179, COMM 180, COMM 181

Use Go Talks for augmentative communication and curriculum enrichment. A powerful, portable, easy to use and very durable speech output device.



### Pocket Go Talk

COMM 82

Pocket Go Talk is a wearable AAC device. With the Pocket Go Talk you can record up to 25 messages, which are easily accessed through five buttons. Each message has 12 seconds of recording time.



### Small Talking Photo Album

COMM 189

With 24 pages and up to 10 seconds of recordable sound for each page, you can use it for daily activity schedules, sequencing, functional reminders and more.



### Large Talking Photo Album

COMM 190

With 38, A4 pages and up to 10 seconds of recordable sound for each page, you can use it for daily activity schedules, sequencing, functional reminders, and more.



### Smart/Speak

COMM 141

AMDI's **Smart/Speak** is a light-weight multi-level augmentative and alternative communication device



### Partner Plus

COMM 157, COM 161



The Partner/Plus is a single message communication device which can also be used alongside other Partner/Plus's. Messages are activated with the slightest movement enabling the user to playback a single recording of up to 30 seconds in length. Can be used with a switch.



### Ipads with Grid Player and Proloquo2go

COMM 145, COMM 163, COMM 165, COMM 166,  
COMM 200 (Mini ipad) COMM 201, COMM 202

Grid Player and Proloquo2go help people who cannot speak or who have unclear speech to communicate. Sentences that are created are spoken out loud. They are suitable for a wide range of people with communication difficulties, including those who use symbols.



### Chatter-Block

COMM 193, COMM 194

A communication aid which is interactive, cross curricular and ideal for speaking and listening. With 6 interactive sides, automatic playback when rolled.



### iPad Rest

COMM 167, COM 168

Clear rest for iPad use.



### Quick Talker

COMM 158, COMM 159

Easy to use go anywhere communicator



### Logan Proxtalker

COMM 129, COMM 137, COMM 138, COMM 199

The Logan Proxtalker is the moveable picture communication device that actually talks and now it can say it all in one go.



### Sequencer

COMM 169

Pressing the button (or a connected switch) causes message one to play. Pressing it again activates message two, then three and so on. This allows individuals to carry on a conversation with a single switch and works well for reading a story, conducting a class report, participating in a school play etc. Sequencer has two message levels for rapid contextual change.



### Big Button with Steps and Levels

COMM 171

A communication aid for single messages or sequences, with multiple levels. 9 levels with 6 seconds recording time plus 1 level with 14 seconds recording time. Two modes – Steps (sequencing) and Levels.



**Cosmonaut**

COMM 210

A large, chunky and grippable stylus for touch screen devices

**AMAneo Bti**

COMM 211

The AMAneo Bti makes it possible to operate an iOS device with any mouse, or assistive mouse including head mouse and thumb mouse.

**Lightwriter SL38**

COMM 191,

Lightwriters are small, robust, portable text-to-speech communication aids specially designed to meet the particular and changing needs of people with speech loss resulting from a wide range of acquired, progressive, and congenital conditions.

**Go Talk Button**

COMM 100, COMM 101

The GoTalk Button records one 10 second message and the large button ensures easy use.



### Lightwriter SL40

COMM 139

In one small package, this lightwriter is an exciting text to speech aid which helps the user to communicate widely and to control their environment.



### Megabee

COMM 116, COMM 117, COMM 118

MegaBee is easy-to-use and rapidly deployed, providing convenient, frequent means of communication with the carer, thus providing a much enhanced lifestyle for the user.



### Partner 4+

COMM 183

A four message communication aid offering vibrating and LED illuminated feedback (optional) for ease of operation by visually impaired users. It combines the high quality of the AMDi Tech Series with portability.



### Big Step by Step

COMM 204

Records a series of messages which the user can step through in sequence. Has interchangeable red, blue, yellow and green tops





### E-Tran Frame & Stand

COMM 212, COMM 213

This frame is portable, sturdy and can be used with people with reasonable eye pointing skills. It comes with removable letters but other symbols can be easily attached to the frame.



### Chatter-Block

COMM 193, COMM 194

A communication aid which is interactive, cross curricular and ideal for speaking and listening. With 6 interactive sides, automatic playback when rolled.



### Partner/One Stepper

COMM 038

The Partner/One Stepper is a single-panel augmentative and alternative communication device with a stepping feature that enables the user to record a sequence of messages for up to a total time of 60 seconds.



### Digital Voice Amplifier

COMM 205

Amplifies low level speech, compact lightweight design, with head worn microphone and throat microphone





## Fast Talk

COMM 206, COMM 207

A pocket sized booklet from the Stroke Association for people with Aphasia to support everyday communication



## I TALK WITH 2 LEVELS

COMM 208

A two buttoned communication device which enables the user to choose between two activities or to ask or answer questions



## Wedge Portable Switch Mount

COMM 209

Enables switches to be placed on an angled surface and the mount has non slip feet so will not move around whilst being used



## TTS Lightweight Choice Board

COMM 020

Wipe able and recordable choice board which enables the user to play recorded messages. Can record up to 4 messages.



**Maxess Mounting Tray**  
**COMM 214**

Enables Switch mounts to be positioned securely to enable comfortable access for the user



**Techcess Allora 2**  
**COMM 215**

The Allora 2 is a keyboard and switch accessible, text to speak communication device which is designed to meet the needs of people who are literate but who have reduced speaking ability.

